Geistlich Fibro-Gide®

Geistlich
Biomaterials

Distributed by:

Geistlich Pharma North America Inc.
Princeton, New Jersey 08540
1-(855) 799-5500
dental.geistlich-na.com

Manufacturer:
Geistlich Pharma AG
CH-6110 Wolhusen
Switzerland

de GERMAN

en ENGLISH

es SPANISH

fr FRENCH

Distributed by:

Geistlich Pharma North America Inc.
Princeton, New Jersey 08540
1-(855) 799-5500
dental.geistlich-na.com

Manufacturer:
Geistlich Pharma AG
CH-6110 Wolhusen
Switzerland

GERMAN

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ENGLISH

COMPOSITION:

Geistlich Fibro-Gide® is a resorbable, porous, non-pyrogenic collagen matrix of porcine origin. The device is of a spongy consistency and is of a white to almost white color. The collagen scaffold is cross linked. The collagen is extracted from veterinary certified pigs and is carefully purified. Geistlich Fibro-Gide® is sterilized in double packaging by Gamma-irradiation. Geistlich Fibro-Gide® is manufactured by a standardized, controlled manufacturing process under an approved quality management system.

PROPERTIES / ACTION:

Geistlich Fibro-Gide® is an implantable device intended for use in soft-tissue augmentation procedures in the oral cavity according to the indications for use. It has mechanical properties appropriate to withstand the mechanical stresses that occur after wound closure in soft-tissue augmentation procedures, i.e. it has good volume stability and it withstands early resorption to allow the formation of new soft-tissue and degrades over time. In addition the matrix is designed with an appropriate thickness to provide sufficient space for the ingrowth of new soft-tissue. Due to its good biological properties, the device becomes well integrated into the surrounding soft-tissue. Geistlich Fibro-Gide® is non-pyrogenic.

Geistlich Fibro-Gide® is provided sterile.

Geistlich Fibro-Gide® is approximately 3 or 6 mm thick. Slight thickness variations with natural collagen materials may be noted.

INDICATIONS FOR USE:

Geistlich Fibro-Gide® is intended for soft-tissue augmentation. It is indicated for:

- > localized gingival augmentation to increase keratinized tissue (KT) around teeth and implants
- > alveolar ridge reconstruction for prosthetic treatment
- > recession defects for root coverage.

CONTRAINDICATIONS:

- > Geistlich Fibro-Gide® must not be used where symptomatic infection or inflammation exists. Before placement, the surgeon should be confident that any symptomatic or recent infection has been properly treated.
- > Geistlich Fibro-Gide® must not be used in patients with a known allergy to collagen.

INSTRUCTIONS FOR USE AND SURGICAL TECHNIQUE:

The general principles of sterile handling and patient treatment must be followed when using Geistlich Fibro-Gide®.

- > The soft-tissue defect located in the oral cavity at the alveolar ridge is exposed. This is achieved by means of (a) flap(s) which allows tension free wound closure.
- > The site is prepared by either creating a split flap or a full thickness flap depending on the site and the available tissue.
- > Geistlich Fibro-Gide® is shaped to the required dimension using a standard sterile surgical scalpel or sterile surgical scissors. To avoid compromising its porous structure, the matrix should not be compressed in dry state.
- > Geistlich Fibro-Gide® is placed on the prepared wound bed under the soft-tissue flap or within the split flap. The placement is most easily achieved with forceps.
- > Take into account that when hydrated the matrix volume will initially expand by approximately 25% corresponding to approximately 5% thickness increase which needs to be taken into account when applying the matrix in a dry state.
- > Complete penetration of the matrix by blood and exudates allows good adaptation and adhesion of the matrix to the surrounding tissues.
- > Once Geistlich Fibro-Gide® is in place, it should be completely covered with a soft-tissue flap that allows healing by primary intention. The soft-tissue flap must be prepared in a way that it can be closed without tension in order to avoid wound dehiscences. Geistlich Fibro-Gide® is usually sufficiently stabilized under the soft-tissue flap. In case additional fixation of Geistlich Fibro-Gide® is required it can be sutured with standard suture techniques using standard sterile suture material and standard sterile suture needles. Caution should be taken when suturing: It is recommended to place the sutures carefully with smooth pulling-forces and not too tight sutures in order to avoid ripping of the matrix.
- > Subsequently primary wound closure is performed by suturing the flap to the surrounding tissue using standard suture techniques.

Geistlich Fibro-Gide® resorbs and is replaced with soft-tissue. Therefore, a second intervention will not be required.

SPECIAL INSTRUCTIONS FOR PERIODONTALLY DISEASED PATIENTS:

A basic requirement for successful periodontal treatment includes eradicating the underlying bacterial infection as well as adequate oral hygiene. Therefore, prior to surgical intervention, patients must receive a hygiene phase of treatment, consisting of oral hygiene instructions, scaling and root planing, and occlusal adjustment when indicated. A postoperative maintenance phase can help to ensure long-term therapeutic success.

POST-OPERATIVE CARE:

The usual post-operative care and medication for soft-tissue augmentation procedures should be followed. Further prosthetic treatment should only be carried out after the healing period, to ensure complete soft-tissue regeneration. Previous studies have indicated that soft-tissue healing, i.e., wound closure and the resolution of normal inflammation, occurs within 4-8 weeks. But clinical judgment, taking into account the healing of both soft and hard tissues in the treated patient, should be applied before prosthetic treatment.

In case of matrix exposure or bacterial contamination, rinsing with appropriate bactericidal solutions is standard practice.

In the rare event where matrix removal is necessary, the tissues adjacent to the matrix should be anesthetized with a local anesthetic. An incision should then be made immediately adjacent to the residual matrix. Following careful reflection of the surrounding tissue, the remaining portion of the matrix can be excised and the area debrided to remove any inflamed or infected tissue.

PRECAUTIONS:

Geistlich Fibro-Gide® must only be used by trained dentists and surgeons.

The matrix should not be used in more than one layer. Geistlich Fibro-Gide® should be used with special caution in patients with:

- > chronic infection at the surgical site (do not use in patients with acute infection at the surgical site)
- > uncontrolled metabolic diseases (e.g. diabetes osteomalacia, thyroid disorder)
- > prolonged corticosteroid therapy
- > autoimmune impairment
- > radiotherapy in the oral cavity
- > healing of soft-tissues might be compromised in patients with insufficiently vascularized soft-tissues (e.g. in smokers).

The material has not been tested in pregnant or lactating women or in children.

ADVERSE REACTIONS:

As the matrix is a collagen product, in very rare cases allergic reactions may occur. In addition the following potential side effects, resulting from the surgical process, may occur: soft-tissue dehiscence, hematoma, increased sensitivity and pain, redness, inflammation and infection.

STORAGE AND HANDLING:

If the sterile packaging is damaged or opened, the product must not be used. The content of the double blister package is designed for single use only.

Do not resterilize Geistlich Fibro-Gide®. Resterilization is not possible. Discard any unused material. Do not use after the expiration date.

Store at controlled room temperature 15°–25 °C (59°–77 °F) and in a dry place. Keep away from sunlight.

The product should be handled using sterile gloves and sterile instruments.

PRESENTATION:

Geistlich Fibro-Gide® is packed in sterile double packaging. Each double packaging contains one matrix.

Re-Order No	Sizes
500838	Geistlich Fibro-Gide® 15 x 20 x 3 mm
500839	Geistlich Fibro-Gide® 20 x 40 x 3 mm
500272	Geistlich Fibro-Gide® 15 x 20 x 6 mm
500273	Geistlich Fibro-Gide® 20 x 40 x 6 mm

Expiry date is printed on the folding box and label of each product.

Distributed by:

Geistlich Pharma North America, Inc.

202 Carnegie Center

Princeton, New Jersey 08540

1 (855) 799-5500

dental.geistlich-na.com

Manufacturer:

Geistlich Pharma AG

Bahnhofstrasse 40

CH-6110 Wolhusen

Switzerland

Made in Switzerland

This product is sold and distributed under US Patent No. 9 655 997 B2

Package insert: dental.geistlich-na.com

CAUTION:

U.S. Federal law restricts this device to sale by or on the order of a dentist or physician.

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