Geistlich Perio-System Combi-Pack

Geistlich
Biomaterials

Geistlich Bio-Gide® Perio

de GERMAN

en ENGLISH

es SPANISH

fr FRENCH

Geistlich Bio-Oss Collagen®

de GERMAN

en ENGLISH

es SPANISH

fr FRENCH

Distributed by:

Geistlich Pharma North America Inc.
Princeton, New Jersey 08540
1-(855) 799-5500
dental.geistlich-na.com

Manufacturer:

Geistlich Pharma AG
CH-6110 Wolhusen
Switzerland

Geistlich Bio-Gide® Perio

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Packungsbeilage: dental.geistlich-na.com

Geistlich Bio-Gide® Perio

Resorbable Bilayer Collagen Membrane for Guided Tissue and Bone Regeneration

Geistlich
Biomaterials

Distributed by:

Geistlich Pharma North America Inc.
Princeton, New Jersey 08540
1-(855) 799-5500
dental.geistlich-na.com

Manufacturer:

Geistlich Pharma AG
CH-6110 Wolhusen
Switzerland

ENGLISH

Geistlich Bio-Gide® Perio is a biocompatible, sterile resorbable bilayer collagen membrane for guided tissue and bone regeneration.

COMPOSITION:

Geistlich Bio-Gide® Perio is a pure collagen membrane obtained by standardized, controlled manufacturing processes. The collagen is extracted from veterinary certified pigs and is carefully purified to avoid antigenic reactions. Geistlich Bio-Gide® Perio is sterilized in double blisters by γ -irradiation. Geistlich Bio-Gide® Perio has a bilayer structure: The porous surface — facing the bone — will allow the ingrowth of the bone forming cells. The dense surface — facing the soft tissue — will prevent the ingrowth of fibrous tissue into the bone defect. The membrane is made of collagen without further cross-linking.

PROPERTIES / ACTION:

The low antigenicity and excellent biocompatibility favor the use of Geistlich Bio-Gide® Perio in dental surgery. Animal studies have shown that the collagen is resorbed in bone cavities within 24 weeks. The fibrous microstructure of Geistlich Bio-Gide® Perio is hydrophilic. The membrane retains its structural integrity even when wet. Fixation by sutures or pins is possible. Geistlich Bio-Gide® Perio has a smooth dense surface. As this prolongs the water uptake, cutting the membrane for periodontal defects is facilitated. The coherent collagen fibers swell and form a unified basic tissue structure. As a result, the adaption to the bone wall and the complete closure of the bony defects is easily achieved. Inflammatory reactions have not been observed, but due to the special physical properties and elongated resorption time a certain inflammatory reaction after placing the membrane cannot be excluded.

INDICATIONS FOR USE:

Geistlich Bio-Gide® Perio is intended for the following uses:

- > augmentation around implants placed in immediate extraction sockets
- > augmentation around implants placed in delayed extraction sockets

- > localized ridge augmentation for later implantation
- > alveolar ridge reconstruction for prosthetic treatment
- > filling of bone defects after root resection, cystectomy, removal of retained teeth
- > guided bone regeneration in dehiscence defects
- > guided tissue regeneration procedures in periodontal defects.

CONTRAINDICATIONS:

Geistlich Bio-Gide® Perio should not be placed where active infection exists. Before placement, the surgeon should be confident that any active or recent infection has been properly treated. Geistlich Bio-Gide® Perio should not be used in patients with a known allergy to collagen.

DIRECTIONS FOR USE:

The general principles of sterile handling and patient medication must be followed when using Geistlich Bio-Gide® Perio:

- > The bone defect is exposed by a mucoperiosteal flap and basic surgical procedures (e.g., curettage) are undertaken.
- > The bone defect is filled with a space-making material, e.g., autogenous bone, Geistlich Bio-Oss®, Geistlich Bio-Oss Collagen®, etc. The defect should not be overfilled.
- > The size of Geistlich Bio-Gide® Perio is adapted by scissors according to the exposed defect. The Geistlich Bio-Gide® Perio package contains sterile stencil templates in 4 different basic periodontal shapes. The membrane should overlap the walls of the defect by at least two millimeters to allow complete bone contact and to prevent gingival connective tissue invasion below the material.
- > The smooth surface should face the soft tissue, and the rough side should face the bone.
- > Geistlich Bio-Gide® Perio is applied over the defect without further treatment and held in place with moderate pressure. The period of time necessary to apply pressure will vary with the degree of bleeding. Adhesion to the bone surface is facilitated by the hydrophilic properties of the membrane.
- > Complete penetration of the membrane by blood and exudate allows perfect adaptation and adhesion of the membrane to bony structures and makes the formation of a blood clot possible.
- > Salivary and other contamination to the material and surgical site should be minimized to avoid bacterial contamination.
- > Due to the high tensile strength of this membrane, fixation is possible. Fixation may be indicated to avoid membrane displacement due to loading or mobilization.
- > The mucoperiosteal flap is sutured over the collagen membrane (e.g., single sutures and deep mattress sutures).
- > Complete wound closure should be attempted to avoid accelerated resorption due to membrane exposure, though clinical experience demonstrates satisfactory healing in most cases when complete wound closure is not possible.
- > In order to avoid the formation of a long junctional epithelium when treating periodontal defects, the coronal portion of Geistlich Bio-Gide® Perio should be adapted closely to the tooth (e.g. with additional fixation using suture material).

POST-OPERATIVE CARE:

- > The patient should be monitored closely.
- > If the membrane becomes exposed, the dehiscence usually heals by itself within several weeks. Membrane removal is usually not necessary. However, to minimize bacterial contamination rinsing with bactericidal solutions is recommended.
- > In the event that early membrane removal is necessary, the tissues adjacent to the membrane should be anesthetized with a local anesthetic. An incision should then be made immediately adjacent to the residual membrane. Following careful reflection of the surrounding tissue, the remaining portion of the membrane can be excised and the area curetted to remove any inflamed or infected tissue.
- > To allow for undisturbed bone regeneration underneath the membrane surgical reentry should not be done before 4 to 6 months postoperatively.

PRECAUTIONS:

Geistlich Bio-Gide® Perio should only be used by trained dentists or oral surgeons. In case of membrane exposure during the healing phase the resorption time may be accelerated. A basic requirement for successful periodontal treatment includes eradicating the underlying bacterial infection as well as adequate oral hygiene. Therefore, prior to surgical intervention, patients must receive a hygiene phase of treatment, consisting of oral hygiene instructions, scaling and root planning, and occlusal adjustment when indicated. A postoperative maintenance phase can help to ensure long-term therapeutic success.

Do not implant the stencil templates. Their perforations enable them to be distinguished from the membrane.

If endosseous implants are involved, the membrane should be used only in combination with a stable implant and not in lieu of achieving primary implant stability. The long-term safety and effectiveness of the membrane in conjunction with bone-filling materials has not yet been established.

Absolute stability of the membrane is important for guided bone regeneration and is a vital condition for therapeutic success, and the smallest movement on the tissue underneath is to be avoided.

Effect on pregnant or lactating women is not known.

Geistlich Bio-Gide® Perio should be used with special caution in patients with

- > acute or chronic infection (e.g., osteomyelitis) at the surgical site
- > uncontrolled metabolic diseases (e.g., diabetes, osteomalacia, thyroid disorder)
- > prolonged corticosteroid therapy
- > autoimmune diseases
- > radiotherapy
- > heavy smoking.

ADVERSE REACTIONS:

As Geistlich Bio-Gide® Perio is a collagen product allergic reactions may not be totally excluded. Possible complications which may occur with any surgery include swelling at the surgical site, flap sloughing, bleeding, dehiscence, hematoma, increased sensitivity and pain, bone loss, redness, and local inflammation.

STORAGE AND HANDLING:

If the sterile packaging is damaged or opened, the product must not be used. The contents of the double blister is designed for single use only.

Do not resterilize Geistlich Bio-Gide® Perio. Resterilization is not possible. Discard any unused material. Do not use after expiry date.

Store at controlled room temperature 15°–25° C (59°–77° F) and in a dry place.

The membrane should be handled using sterile gloves or sterile instruments.

PRESENTATION: Geistlich Bio-Gide® Perio is packed in sterile double blisters.

Re-order #	Size
20148	1 Geistlich Bio-Gide® Perio membrane, 16 x 22 mm

Geistlich Perio-System Combi-Pack

Re-order #	Size
20147	1 Geistlich Bio-Gide® Perio membrane, 16 x 22 mm 1 Geistlich Bio-Oss Collagen®, 100 mg (1 block 100 mg Geistlich Bio-Oss® spongiosa (cancellous) granules and 10% collagen

Distributed by:

Geistlich Pharma North America Inc.
Princeton, New Jersey 08540
1-(855) 799-5500
dental.geistlich-na.com

Manufacturer:

Geistlich Pharma AG
CH-6110 Wolhusen
Switzerland

Made in Switzerland

CAUTION: Federal law restricts this device to sale by or on the order of a dentist or physician.
Package insert: dental.geistlich-na.com

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Geistlich Pharma AG
CH-6110 Wolhusen
Switzerland

ENGLISH

Geistlich Bio-Oss Collagen® is a sterile, biocompatible bone mineral plus collagen for bone grafting in periodontal, oral and maxillofacial surgery.

COMPOSITION:

Geistlich Bio-Oss Collagen® is a mixture of spongiosa (cancellous) bone granulate and 10% porcine collagen fibers in a block form, and is sterilized by γ -irradiation.

PROPERTIES / ACTION:

The anorganic bone matrix of Geistlich Bio-Oss Collagen® has physical and chemical properties comparable to the mineralized matrix of human bone. The collagen facilitates handling of the graft particles and acts to hold the Geistlich Bio-Oss® at the desired place. The consistency of this material readily allows it to take the shape of the defect. The graft is slowly resorbed and replaced by new bone.

INDICATIONS FOR USE:

Geistlich Bio-Oss Collagen® is intended for the following uses:

- > augmentation or reconstructive treatment of the alveolar ridge
- > filling of periodontal defects
- > filling of defects after root resection, apicoectomy, and cystectomy
- > filling of extraction sockets to enhance preservation of the alveolar ridge
- > elevation of the maxillary sinus floor
- > filling of periodontal defects in conjunction with products intended for Guided Tissue Regeneration (GTR) and Guided Bone Regeneration (GBR)
- > filling of peri-implant defects in conjunction with products intended for Guided Bone Regeneration (GBR).

CONTRAINDICATIONS:

Geistlich Bio-Oss Collagen® should not be used in the presence of infected wounds. Geistlich Bio-Oss Collagen® should not be used in patients with a known allergy to collagen.

DIRECTIONS FOR USE:

The general principles of sterile handling and patient medication must be followed when using Geistlich Bio-Oss Collagen®.

- > The filling of periodontal defects with Geistlich Bio-Oss Collagen® requires successful treatment of the periodontal lesion (e.g., root planing, debridement of granulation tissue) prior to implantation. The defect should be covered with a membrane (e.g., Geistlich Bio-Gide®) for optimal tissue regeneration.
- > After exposure of the bony defect with mucoperiosteal flap, all granulation tissue must be carefully removed.
- > The material should be placed in direct contact with well vascularized, bleeding bone surfaces. Cortical bone should be mechanically perforated to facilitate ingrowth of new blood vessels and bone forming cells.
- > Application: Geistlich Bio-Oss Collagen® may be cut to the required size with sterile tweezers and scissors, in dry form or after moistening with blood or physiological saline solution. The material is applied using sterile instruments (e.g., tweezers).
- > If necessary, the material can be slightly molded in situ using a spatula or similar instrument.
- > Overfilling of the defect should be avoided.
- > It is advisable to cover Geistlich Bio-Oss Collagen® with a membrane barrier (e.g., Geistlich Bio-Gide®).
- > When closing the wound, the soft tissue flap should completely cover the implanted Geistlich Bio-Oss Collagen®, and should be fixed by sutures.
- > If primary wound closure cannot be achieved completely, further mobilization of the flap (incision through the periosteum) should be performed.
- > A surgical dressing may be placed over the wound for one to two weeks.
- > For improved bone formation in large osseous defects with one or two bony walls, Geistlich Bio-Oss Collagen® may be mixed with autogenous particulate spongiosa (cancellous) bone and marrow or covered with a cell occlusive membrane (e.g., Geistlich Bio-Gide®).

PRECAUTIONS:

Geistlich Bio-Oss Collagen® should only be used by trained dentists or oral surgeons.

In order to facilitate the formation of new bone, Geistlich Bio-Oss Collagen® should only be implanted in direct contact with well vascularized bony tissue (selective osteoplasty of adjacent cortical bone may be necessary) and if necessary microfracturation of the adjacent bone might be performed.

From experience, mechanical loading (compression loading) of Geistlich Bio-Oss augmented areas is possible after 6 months at the earliest. Importantly the appropriate timing of implant insertion depends on the local bone volume at the time of implant insertion.

In larger defects a mixture of autogenous bone or bone marrow may improve the formation of new bone.

Geistlich Bio-Oss Collagen® should be used with special caution in patients with:

- > acute or chronic infection (osteomyelitis) at the surgical site
- > metabolic diseases (diabetes, hyperparathyroidism, osteomalacia)
- > osteoporosis
- > severe renal dysfunction, severe liver disease
- > high dose corticosteroid therapy
- > vascular impairment at the implant site
- > autoimmune disease
- > radiotherapy
- > heavy smoking.

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ADVERSE REACTIONS:

Incompatibility reactions with Geistlich Bio-Oss Collagen® cannot be excluded. As Geistlich Bio-Oss Collagen® contains collagen, in very rare cases allergic reactions may occur. Possible complications which may occur with any surgery include swelling at the surgical site, flap sloughing, bleeding, local inflammation, bone loss, infection or pain.

STORAGE AND HANDLING:

If the sterile packaging is damaged or opened the product must not be used. The contents of the blister package are designed for single use only. Discard any unused material after opening.

Resterilization is not possible.

Do not use after the expiration date.

Store at controlled room temperature 15°–25 °C (59°–77 °F) and in a dry place.

PRESENTATION: Geistlich Bio-Oss Collagen®:

Re-Order #:	Packaged:
500413	1 block 50 mg Geistlich Bio-Oss® spongiosa (cancellous) granules and 10% collagen
20141	1 block 100 mg Geistlich Bio-Oss® spongiosa (cancellous) granules and 10% collagen
20142	1 block 250 mg Geistlich Bio-Oss® spongiosa (cancellous) granules and 10% collagen
20143	1 block 500 mg Geistlich Bio-Oss® spongiosa (cancellous) granules and 10% collagen

Geistlich Combi-Kit Collagen®

Re-Order #:	Packaged:
20146	1 block 100 mg Geistlich Bio-Oss® spongiosa (cancellous) granules and 10% collagen, and 1 Geistlich Bio-Gide® membrane, 16 x 22 mm

Geistlich Perio-System Combi-Pack

Re-Order #:	Packaged:
20147	1 block 100 mg Geistlich Bio-Oss® spongiosa (cancellous) granules and 10% collagen, and 1 Geistlich Bio-Gide® Perio membrane, 16 x 22 mm

Distributed by:

Geistlich Pharma North America Inc.

Princeton, New Jersey 08540

1-(855) 799-5500

dental.geistlich-na.com

Manufacturer:

Geistlich Pharma AG

CH-6110 Wolhusen

Switzerland

Made in Switzerland

CAUTION: Federal law restricts this device to sale by or on the order of a dentist or physician.

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